



Khaw Surgical Gonioprism

	Product Code	Contact OD	<i>Designed with: P. T Khaw, M.D. PhD London, U.K.</i>	
	OKSG CE	11.5mm		

Lens Design

- § The Khaw surgical gonioprism was designed for easy manipulation during goniotomy and direct viewing gonioscopy procedures.
- § A fixation ring attached to the posterior portion of the lens helps to stabilize the globe during surgical procedures.
- § Wide field of view lens provides a clear view of the anterior chamber and anterior chamber angle during goniotomy surgical procedures.
- § Black anodized aluminum handle provides easy manipulation of the lens.
- § Lens magnification provides detailed view of the anterior chamber.
- § Use of a coupling viscoelastic may be of considerable help visualizing the anterior chamber angle and to prevent air bubbles in the lens corneal interface.

Cleaning and Disinfection

See Cleaning Method 1

