





Ocular Ahmed 1.5X Surgical Gonio Lens

	Product Code	Gonio Mag	Static FOV	Contact OD	Ring Diameter	Handle Length	
	OASG	1.50x	90°	9.90mm	28.6mm	-	
	OASG-H	1.50x	90°	9.90mm	26.0mm	72mm	
							

Lens Design

§ Designed to allow viewing of anterior chamber and chamber angle during anterior chamber angle surgical procedures

§ Excellent for quick angle visualization with no need for patient or microscope tilt.

§ Can also be used for procedures like goniosynechialysis, and certain suprachoroidal and Schlemm's canal device implantations and other ab-interno devices.

Cleaning and Sterilization

See Cleaning Method 3



2255 116th Ave NE, Bellevue, Washington 98004-3039 USA
 T: 425-455-5200 or 800-888-6616 F: 425-462-6669
 E: ocular@ocularinc.com I: www.ocularinc.com

© 2012 Ocular Instruments
 11027-3483